

EZ- Nitric Oxide Assay Kit

Metabolism assay kit
(Colorimetric)

Cat. No. DG-NO500

FOR RESEARCH USE ONLY.

NOT FOR USE IN DIAGNOSTIC PROCEDURES.

▪ Product Description

Nitric oxide (NO) is a radical compound generated during the conversion of L-arginine to citrulline by nitric oxide synthase (NOS). It plays a critical role in neurotransmission, blood vessel regulation, immune responses, and cell death. NO also influences macrophage-mediated cytotoxicity and is known to inhibit mitochondrial respiration and DNA synthesis.

The EZ-Nitric Oxide Assay Kit is based on the diazotization reaction (Griess method) and measures the NO concentration in the sample. This kit accurately quantifies nitrate (NO_3^-) levels by measuring nitrite (NO_2^-), a conversion product of nitrate present in tissues or cells. The reaction between sulfanilamide and nitrite, followed by the addition of *N*-(1-naphthyl)ethylenediamine, generates a product that can be used to quantify the total NO_3^- concentration.




Fig 1. Nitrite quantification using the Griess reagent.

▪ Contents and Storage Conditions

Component	500 assay	Storage
Griess Reagent A	25 mL	4 °C
Griess Reagent B	25 mL	
Nitrite Standard (2 mM)	500 µL X 2 vial	

* It remains stable for 6 months at 4°C after opening.

* Do not mix Griess reagent A and B in advance; mix them just before use.

▪ General Protocol

Sample preparation

Note

- If the sample contains proteins, filter it using a 0.2 μm filter before using it in the experiment.

● **Serum/Plasma**

- ① Use a filter to remove proteins from the sample.
- ② Centrifuge the sample at $10,000 \times g$ for 10 minutes at 4°C and use only the supernatant for analysis.

● **Urine**

- ① Dilute the sample at least 10 times using PBS, filter it, and use it immediately.
- ② For more accurate measurements, prepare samples at various concentrations using PBS.

● **Cell lysate**

- ① Homogenize 1×10^6 cells quickly with 100 μl of ice-cold PBS.
- ② Incubate on ice for 10 minutes, then centrifuge at $10,000 \times g$ for 5 minutes at 4°C. Transfer the supernatant to a new tube.
- ③ Add 10–100 μl of the sample to each well, then add PBS to adjust the total volume to 100 μl for use in the experiment.

● **Culture supernatant**

- ① Avoid using media that contain nitrite.
- ② Transfer 100 μl of the culture supernatant (without filtration) to a 96-well plate, and use it directly as the sample for analysis.
- ③ When using media with a strong color, dilute the nitrite standard solution in the same medium and prepare the standard curve.

Nitrite standard preparation

: Prepare the nitrite standard solution by diluting 2 mM nitrite solution and D.W as follows.

Standard Tube No.	Nitrite standard solution	D.W	Final standard Nitrite Conc. (μM)
1	50 μl of 2 mM (2000 μM)	450 μl	200
2	250 μl of Tube # 1	250 μl	100
3	250 μl of Tube # 2	250 μl	50
4	250 μl of Tube # 3	250 μl	25
5	250 μl of Tube # 4	250 μl	12.5
6	250 μl of Tube # 5	250 μl	6.25
7	250 μl of Tube # 6	250 μl	3.125
8	0 μl	250 μl	0

* It is recommended that the standard be measured for each experiment.

* Prepare the Nitrite standard solution and sample in a volume of at least 250 μl to allow for testing in duplicate or more for each concentration.

Nitric Oxide (Griess Reagent) Assay

- ① Before the experiment, store Griess reagents A and B at room temperature for about 15–30 minutes before use.
- ② Add 100 μl of the Nitrite standard solution and sample into each well according to the concentration.
- ③ Add 50 μl of Griess reagent A to each well and incubate in the dark at room temperature for 10 minutes.
- ④ Add 50 μl of Griess reagent B to each well and incubate again in the dark at room temperature for 10 minutes.
- ⑤ After the reaction, measure the absorbance at 540 nm using a microplate reader. If a 540 nm wavelength is unavailable, measure within the range of 520–550 nm.

▪ Calculation

Nitrite concentration of sample (μM)

1. Correct the O.D. values for all measurements by subtracting the blank's O.D. value.
* Blank = O.D._{540nm} from Nitrite standard No. 8 (0 μM Nitrite)
2. Plot the standard curve using the O.D. values (Y-axis) versus the concentrations of Nitrite standard (X-axis) to obtain the linear equation.
3. Calculate the Nitrite amount in the sample by substituting the sample's O.D. value into the standard curve.
* If a background control has been set, subtract the background control measurement from the sample's measurement before calculating the Nitrite amount.
* If the sample was diluted 2-fold, multiply the calculated Nitrite concentration (μM) by 2 to obtain the accurate Nitrite amount.




Fig 2. Nitrite standard curve. Assay was performed following the kit protocol.

▪ Related Product

	Products	Catalog No.	Assay
Oxidative Stress Assay Kit	EZ-Superoxide Dismutase (SOD) Assay Kit (Colorimetric)	DG-SOD400	400 Assay
	EZ-Glutathione Assay Kit (Colorimetric)	DG-GLU200	200 Assay
	EZ-Catalase Assay Kit (Fluorometric/Colorimetric)	DG-CAT400	400 Assay
	EZ-Hydrogen peroxide/Peroxidase Assay Kit (Fluorometric/Colorimetric)	DG-PER500	500 Assay
	EZ-Lipid Peroxidation (TBARS) Assay Kit (Colorimetric)	DG-TBA200	200 Assay
	EZ-Total Antioxidant Capacity (TAC) Assay Kit (Colorimetric)	DG-TAC200	200 Assay
	EZ-DPPH Antioxidant Assay Kit (Colorimetric)	DG-DPH400	400 Assay
Metabolism Assay Kit	EZ-Lactate Assay Kit (Colorimetric)	DG-LAC100	100 Assay
	EZ-Acetylcholinesterase Assay Kit (Colorimetric)	DG-ACE100	100 Assay
	EZ-Ascorbic Acid Assay Kit (Colorimetric)	DG-ASC100	100 Assay
	EZ-ATP Assay Kit (Fluorometric/Colorimetric)	DG-ATP100	100 Assay
	EZ-Free Fatty Acid Assay Kit (Fluorometric/Colorimetric)	DG-FFA100	100 Assay
	EZ-Free Glycerol Assay Kit (Fluorometric/Colorimetric)	DG-FGC100	100 Assay
	EZ-Glucose Assay Kit (Fluorometric/Colorimetric)	DG-GCS100	100 Assay
	EZ-HDL, LDL/VLDL Assay Kit (Fluorometric/Colorimetric)	DG-CHO100	100 Assay
	EZ-Total Cholesterol Assay Kit (Fluorometric/Colorimetric)	DG-TSC100	100 Assay
	EZ-Triglyceride Quantification Assay Kit (Fluorometric/Colorimetric)	DG-TGC100	100 Assay